

Perfusion Fluid CNS Dextran for Microdialysis

Intended Purpose

The Perfusion Fluid CNS Dextran is intended to act as a carrier inside the CMA Microdialysis Probes that have 55 kDa and 100 kDa membranes.

Ordering Information

	Ref. No.	Package Size
Perfusion Fluid CNS Dextran	8050151	Glass bottles, 2 x 7.5 mL



Technical Information

Features	Isotonic, Sterile
Contents	NaCl 147 mmol/L KCL 2.7 mmol/L CaCl ₂ 1.2 mmol/L MgCl ₂ 0.85 mmol/L Total Chloride Contents: 153.8 mmol/L Dextran MW 200 kDa 3%
pH	~6

Symbol Explanation

	Fulfils EU Medical Device Regulation (MDR) 2017/745
	Sterilized by use of aseptic processing technique
	Single use only
	Storage temperature: 15 - 25 °C
	Medical Device