


Perfusion Fluid CNS for Microdialysis

Intended Purpose

The Perfusion Fluid CNS is intended to act as a carrier inside the CMA Microdialysis Probes used in the brain.

Ordering Information

	Ref. No.	Package Size
Perfusion Fluid CNS	P000151	Glass bottles, 10 x 7.5 mL



Technical Information

Features	Isotonic, Sterile
Contents	NaCl 147 mmol/L KCL 2.7 mmol/L CaCl ₂ 1.2 mmol/L MgCl ₂ 0.85 mmol/L Total Chloride Contents: 153.8 mmol/L
Osmolality	305 mOsm/kg
pH	~6

Symbol Explanation

	Fulfills EU Medical Device Regulation (MDR) 2017/745
	Sterilized using steam
	Single use only
	Storage temperature: 15 - 25 °C
	Medical Device